

Pegylated Remote Loadable Liposome (2mL vial)

FOR RESEARCH USE ONLY. NOT INTENDED FOR HUMAN THERAPEUTIC OR DIAGNOSTIC USE

300202
Page 1 of 1

Manufactured by:	Avanti Polar Lipids, Inc.
Formulation:	A pegylated, long-circulating liposomal formulation with an ammonium sulfate core for remote loading. Liposomes are stored frozen in a 10% w/v sucrose 10mM histidine buffer.
Physical form:	Frozen translucent liposomal dispersion
Storage:	Freezer (-20° C)
Expiration date:	One year from date of receipt
M Lot Number:	6540SAA010
Lot Number:	300202S-1EA-A-010

ANALYSIS	SPECIFICATION	RESULTS
Total Lipid Content (HPLC)	Run & Report	25.0 mg/mL
Lipids Ratio (HPLC):		
Cholesterol	Run & Report	4.9 mg/mL
DSPE-mPEG 2000	Run & Report	4.7 mg/mL
HSPC	Run & Report	15.4 mg/mL
Particle size (DLS)	Run & Report	88.4 nm

Approved By: Signature On File
 Reviewed By: Signature On File