

**Pegylated Remote Loadable Liposome (5mL vial)**
**300205S-1ea**
**Page 1 of 1**

<b>Physical State</b>	Liposomal dispersion
<b>Storage</b>	-20°C
<b>Expiration Date</b>	One year from date of receipt
<b>M Lot Number</b>	6541SAA010
<b>Avanti Lot Number</b>	300205S-1EA-010

ANALYSIS	SPECIFICATION	RESULTS
Physical Examination	Translucent liposomal dispersion	Pass
HPLC-ELSD (total lipid content)	Run & Report (mg/mL)	25.0 mg/mL
HPLC-ELSD (lipids ratio)		
<i>Cholesterol</i>	Run & Report (%)	4.9 mg/mL (19.6%)
<i>DSPE-mPEG 2000</i>	Run & Report (%)	4.7 mg/mL (18.8%)
<i>HSPC</i>	Run & Report (%)	15.4 mg/mL (61.6%)
Particle Size (DLS)	Run & Report	88.4 nm

Approved By: *Melissa Liles*

For Research Use Only. Not Intended For Human Therapeutic Or Diagnostic Use.