

GD3 ELISA (Clear)

 GD3 ELISA Snoopers[®], Clear, for use in absorbance applications

330657-1ea

Page 1 of 1

Storage	-20°C
Expiration Date	One year from date of receipt
M Lot Number	5039NAA011
Avanti Lot Number	330657-1EA-011

ANALYSIS	SPECIFICATION	RESULTS
Confirm release of GD3 (Bovine Milk)	Physical Examination Powder: White solid which contains no foreign matter	Pass
	TLC (5:4:1 Chloroform: Methanol: Water w/ MgCl ₂) >99% Purity Ninhydrin: negative Iodine: one major spot Phosphorus: negative Charring: positive Water dip: one major spot	All Pass
	Proton NMR NMR spectrum consistent with structure	Consistent with structure
	Mass Spectroscopy Consistent with structure	Consistent with structure
Lipid Content	100µL of 10µM GD3 solution applied to each well and evaporated using proprietary technique	Completed
Packaging	Packaged in foil bag under argon gas and heat sealed	Completed

Approved By:



For Research Use Only. Not For Human Use.