

## Soy Phospholipid Mixture

Phospholipids from Soybean, Contains Soy PC (3.8mg) Soy PE (3.0mg) Soy PI (1.8mg) Soy PA (0.7mg) Soy LPC (0.7mg) with 0.1% BHT

690050C-10MG

Page 1 of 1

<b>Physical State</b>	Chloroform
<b>Storage</b>	-20°C
<b>Expiration Date</b>	One year from date of receipt
<b>M Lot Number</b>	5083CIA017
<b>Avanti Lot Number</b>	690050C-10MG-A-017

ANALYSIS	SPECIFICATION	RESULTS
Physical examination	<b>Chloroform:</b> Clear solution with no foreign matter	Pass
Quantitative Phosphorus NMR	Run and report composition per 10mg sample: 38:30:18:7:7 (PC:PE:PI:PA:LPC)	37.4 : 30.9 : 17.7 : 7.1 : 6.9 % PC : PE : PI : PA : LPC

Approved By:

