

Soy Phospholipid Mixture

Phospholipids from Soybean, Contains Soy PC (3.8mg) Soy PE (3.0mg) Soy PI (1.8mg) Soy PA (0.7mg) Soy LPC (0.7mg) with 0.1% BHT

690050P-10MG

Page 1 of 1

Physical State	Powder
Storage	-20 °C
Expiration Date	One year from date of receipt
M Lot Number	5083PIA016
Avanti Lot Number	690050P-10MG-A-016

ANALYSIS	SPECIFICATION	RESULTS
Physical examination	Powder: White solid which contains no foreign matter	Pass
Quantitative Phosphorus NMR	Run and report composition per 10mg sample: 38:30:18:7:7 (PC:PE:PI:PA:LPC)	37.4:31.1:18.0:7.0:6.5 % PC:PE:PI:PA:LPC

Approved By:

