

Soy Phospholipid Mixture

Phospholipids from Soybean, Contains Soy PC (3.8 mg), Soy PE (3.0 mg), Soy PI (1.8 mg), Soy PA (0.7 mg), Soy LPC (0.7 mg) with 0.1% BHT

690050C-10mg

Page 1 of 1

| | |
|-------------------|-------------------------------|
| Physical State | Chloroform |
| Storage | -20°C |
| Expiration Date | One year from date of receipt |
| M Lot Number | 5083CIA018 |
| Avanti Lot Number | 690050C-10MG-A-018 |

| ANALYSIS | SPECIFICATION | RESULTS |
|-----------------------------|---|--|
| Physical Examination | Chloroform: Clear solution with no foreign matter | Pass |
| Quantitative Phosphorus NMR | Run and report composition per 10mg sample: 38:30:18:7:7 (PC: PE: PI: PA: LPC) | 38.43% PC 30.70% PE 17.22% PI 7.04% PA 6.61% LPC |

Approved By:



For Research Use Only. Not For Human Use.