

Soy Phospholipid Mixture

Phospholipids from Soybean, Contains Soy PC (3.8mg), Soy PE (3.0mg), Soy PI (1.8mg), Soy PA (0.7mg), Soy LPC (0.7mg) with 0.1% BHT

690050C-10mg**Page 1 of 1**

Physical State Chloroform
Storage -20°C
Expiration Date One year from date of receipt
M Lot Number 5083CIA021
Avanti Lot Number 690050C-10MG-A-021

ANALYSIS	SPECIFICATION	RESULTS
Physical Examination	Chloroform: Clear solution which contains no foreign matter	Pass
Quantitative Phosphorus NMR	Run and report composition per 10mg sample: 38:30:18:7:7 (PC:PE:PI:PA:LPC)	3.93 mg PC 3.12 mg PE 1.85 mg PI 0.61 mg PA 0.79 mg LPC

Approved By:



For Research Use Only. Not For Human Use.