

Soy Phospholipid Mixture

Phospholipids from Soybean, Contains Soy PC (3.8mg), Soy PE (3.0mg), Soy PI (1.8mg), Soy PA (0.7mg), Soy LPC (0.7mg) with 0.1% BHT

690050C-10MG

Page 1 of 1

Physical State	Chloroform
Storage	-20°C
Expiration Date	One year from date of receipt
M Lot Number	5083CIB017
Avanti Lot Number	690050C-10MG-B-017

ANALYSIS	SPECIFICATION	RESULTS
Physical Examination	Chloroform: Clear solution with no foreign matter	Pass
Quantitative Phosphorus NMR	Run and report composition per 10mg sample: 38:30:18:7:7 (PC:PE:PI:PA:LPC)	37.4:30.9:17.7:7.1:6.9% PC:PE:PI:PA:LPC

Approved By:



For Research Use Only. Not For Human Use.