

Soy Phospholipid Mixture

Phospholipids from Soybean, Contains Soy PC (3.8 mg), Soy PE (3.0 mg), Soy PI (1.8 mg), Soy PA (0.7 mg), Soy LPC (0.7 mg) with 0.1% BHT

690050C-10mg

Page 1 of 1

Physical State Chloroform
Storage -20°C
Expiration Date One year from date of receipt
M Lot Number 5083CIB018
Avanti Lot Number 690050C-10MG-B-018

ANALYSIS	SPECIFICATION	RESULTS
Physical Examination	Chloroform: Clear solution with no foreign matter	Pass
Quantitative Phosphorus NMR	Run and report composition per 10mg sample: 38:30:18:7:7 (PC: PE: PI: PA: LPC)	38.43% PC 30.70% PE 17.22% PI 7.04% PA 6.61% LPC

Approved By:



For Research Use Only. Not For Human Use.