

Chol-PEG600

Cholesterol-(polyethylene glycol-600)

880001M-25mg

Page 1 of 1

Molecular Weight	1003.389 (average MW due to polydispersity of PEG)
Physical State	Methylene Chloride
Storage	-20 °C
Expiration Date	One year from date of receipt
M Lot Number	6282MJA010
Avanti Lot Number	880001M-25MG-A-010

ANALYSIS	SPECIFICATION	RESULTS
Physical examination	Methylene Chloride: Clear solution with no foreign matter.	Pass
TLC (9:1 Chloroform: Methanol)	>99% Purity Ninhydrin: negative Iodine: one major spot Phosphorus: negative Charring: positive Water dip: one major spot	All Pass
Proton NMR	NMR spectrum consistent with structure	Consistent with structure

Approved By:

