

DSPE-PEG(2000) Cyanur

1,2-distearoyl-sn-glycero-3-phosphoethanolamine-N-[cyanur(polyethylene glycol)-2000] (ammonium salt)

880122P-10mg

Page 1 of 1

Molecular Weight	2938.436 (average MW due to polydispersity of PEG)
Physical State	Powder
Storage	-20 °C
Expiration Date	One year from date of receipt
M Lot Number	6285PIA015
Lot Number	880122P-10MG-A-015

ANALYSIS	SPECIFICATION	RESULTS
Physical examination	Powder: White solid which contains no foreign matter.	Pass
TLC (85:15 Chloroform: Methanol)	>99% Purity Ninhydrin: negative Iodine: one major spot Phosphorus: positive Charring: negative Water dip: one major spot	All Pass
Proton NMR	NMR spectrum consistent with structure	Consistent with structure

Approved by:

