

18:1 PEG2000 PE

1,2-dioleoyl-sn-glycero-3-phosphoethanolamine-N-[methoxy(polyethylene glycol)-2000] (ammonium salt)

880130P-1g
Page 1 of 1

| | |
|-------------------------|--|
| Molecular Weight | 2801.465 (average MW due to polydispersity of PEG) |
| Physical State | Powder |
| Storage | -20 °C |
| Expiration Date | One year from date of receipt |
| M Lot Number | 6292PQA047 |
| Lot Number | 880130P-1G-A-047 |

| ANALYSIS | SPECIFICATION | RESULTS |
|-------------------------------------|---|-----------------|
| Physical examination | Powder: White solid which contains no foreign matter. | Pass |
| TLC (85:15 Chloroform: Methanol) | >99% Purity Ninhydrin: negative Iodine: one major spot Phosphorus: positive Charring: positive Water dip: one major spot | All pass |
| HPLC-ELSD | NLT 99% (AUC) | 99% (AUC) |
| FAME (GC/FID) | NLT 99% 18:1 (AUC) | 100% 18:1 (AUC) |

Approved by:

