

DSPE-PEG(2000) Carboxylic Acid

1,2-distearoyl-sn-glycero-3-phosphoethanolamine-N-[carboxy(polyethylene glycol)-2000] (sodium salt)

880135P-100mg

Page 1 of 1

Molecular Weight	2780.38 (average MW due to polydispersity of PEG)
Chemical Formula	C ₁₃₁ H ₂₅₇ NNaO ₅₅ P (average formula due to polydispersity of PEG)
Physical State	Powder
Storage	-20°C
Expiration Date	One year from date of receipt
M Lot Number	6296PLD015
Lot Number	880135P-100MG-D-015

ANALYSIS	SPECIFICATION	RESULTS
Physical Examination	Powder: White solid which contains no foreign matter	Pass
TLC (80:20:2 Chloroform: Methanol: Ammonium Hydroxide)	NLT 98% Purity Ninhydrin: negative Iodine: one major spot Phosphorus: positive Charring: negative Water dip: one major spot	All Pass
HPLC-ELSD	NLT 98% (AUC)	98% DSPEPEG2KCOOH (AUC) 2% MPEG2K
FAME (GC/FID)	NLT 99% 18:0 (AUC)	100% (AUC)
Proton NMR	NMR spectrum consistent with structure	Consistent with structure
Molecular Weight (Mw) by QTOF-MS	Run and Report	3070.720
Polydispersity Index (PDI) by QTOF-MS	Run and Report	1.002

 Approved by: *Anderson Thomas*

For Research Use Only. Not For Human Use.