

## DSPE-PEG(2000) Carboxylic Acid

1,2-distearoyl-sn-glycero-3-phosphoethanolamine-N-[carboxy(polyethylene glycol)-2000] (sodium salt)

880135P-50mg

Page 1 of 1

<b>Molecular Weight</b>	2780.38 (average MW due to polydispersity of PEG)
<b>Chemical Formula</b>	C <sub>131</sub> H <sub>257</sub> NNaO <sub>55</sub> P (average formula due to polydispersity of PEG)
<b>Physical State</b>	Powder
<b>Storage</b>	-20°C
<b>Expiration Date</b>	One year from date of receipt
<b>M Lot Number</b>	6296PKB011
<b>Lot Number</b>	880135P-50MG-B-011

ANALYSIS	SPECIFICATION	RESULTS
Physical examination	Powder: White solid which contains no foreign matter.	Pass
TLC (80:20:2 Chloroform:Methanol:Ammonium Hydroxide)	>99% Purity Ninhydrin: negative Iodine: one major spot Phosphorus: positive Charring: negative Water dip: one major spot	All Pass
HPLC-ELSD	NLT 99% (AUC)	97.0% (AUC) 3.0% MPEG2000 (AUC)
FAME (GC/FID)	NLT 99% 18:0 (AUC)	100% (AUC)
Proton NMR	NMR spectrum consistent with structure	Consistent with structure

Approved by:

