

**18:0 PEG5000 PE**

1,2-distearoyl-sn-glycero-3-phosphoethanolamine-N-[methoxy(polyethylene glycol)-5000] (ammonium salt)

**880220P-200mg**
**Page 1 of 1**

<b>Molecular Weight</b>	5801.071 (average MW due to polydispersity of PEG)
<b>Physical State</b>	Powder
<b>Storage</b>	-20°C
<b>Expiration Date</b>	One year from date of receipt
<b>M Lot Number</b>	6303PNA050
<b>Lot Number</b>	880220P-200MG-A-050

ANALYSIS	SPECIFICATION	RESULTS
Physical examination	<b>Powder:</b> White solid which contains no foreign matter.	Pass
TLC	>99% Purity Ninhydrin: negative Iodine: one major spot Phosphorus: positive Charring: negative Water dip: one major spot	All Pass
HPLC-ELSD	NLT 99% (AUC)	99.9% PEG2PE 0.09% MPEG2000
FAME (GC/FID)	NLT 99% 18:0 (AUC)	100.0% 18:0 (AUC)

Approved by:

