

18:0 PEG5000 PE

1,2-distearoyl-sn-glycero-3-phosphoethanolamine-N-[methoxy(polyethylene glycol)-5000] (ammonium salt)

880220P-500mg

Page 1 of 1

Molecular Weight	5801.071 (average MW due to polydispersity of PEG)
Physical State	Powder
Storage	-20°C
Expiration Date	One year from date of receipt
M Lot Number	6303PPA049
Lot Number	880220P-500MG-A-049

ANALYSIS	SPECIFICATION	RESULTS
Physical examination	Powder: White solid which contains no foreign matter.	Pass
TLC	>99% Purity Ninhydrin: negative Iodine: one major spot Phosphorus: positive Charring: negative Water dip: one major spot	All Pass
HPLC-ELSD	NLT 99% (AUC)	99.5% (AUC)
FAME (GC/FID)	NLT 99% 18:0 (AUC)	100.0% 18:0 (AUC)

Approved by:

