

18:1 PEG5000 PE

1,2-dioleoyl-sn-glycero-3-phosphoethanolamine-N-[methoxy(polyethylene glycol)-5000] (ammonium salt)

880230C-25mg
Page 1 of 1

Molecular Weight	5797.040 (average MW due to polydispersity of PEG)
Physical State	Chloroform
Storage	-20°C
Expiration Date	One year from date of receipt
M Lot Number	6309CJD034
Avanti Lot Number	880230C-25MG-D-034

ANALYSIS	SPECIFICATION	RESULTS
Physical Examination	Chloroform: Clear solution with no foreign matter	Pass
TLC (85:15 Chloroform: Methanol)	>99% Purity Ninhydrin: negative Iodine: one major spot Phosphorus: positive Charring: positive Water dip: one major spot	All Pass
HPLC-ELSD	NLT 99% (AUC)	99.5% (AUC)
FAME (GC/FID)	NLT 99% 18:1 (AUC)	100% (AUC)
Proton NMR	NMR spectrum consistent with structure	Consistent with structure

 Approved By: 

For Research Use Only. Not For Human Use.