

16:0 PEG3000 PE

1,2-dipalmitoyl-sn-glycero-3-phosphoethanolamine-N-[methoxy(polyethylene glycol)-3000] (ammonium salt)

880300C-25mg
Page 1 of 1

Molecular Weight	3718.547 (average MW due to polydispersity of PEG)
Physical State	Chloroform
Storage	-20°C
Expiration Date	One year from date of receipt
M Lot Number	6312CJB012
Lot Number	880300C-25MG-B-012

ANALYSIS	SPECIFICATION	RESULTS
Physical examination	Chloroform: Clear solution with no foreign matter.	Pass
TLC (85:15 Chloroform: Methanol)	>99% Purity Ninhydrin: negative Iodine: one major spot Phosphorus: positive Charring: negative Water dip: one major spot	All Pass
HPLC-ELSD	NLT 99% (AUC)	99.6% (AUC)
FAME (GC/FID)	NLT 99% 16:0 (AUC)	100% (AUC)

Approved by:

