

14:0 PEG3000 PE

1,2-dimyristoyl-sn-glycero-3-phosphoethanolamine-N-[methoxy(polyethylene glycol)-3000] (ammonium salt)

880310P-25mg

Page 1 of 1

Molecular Weight	3662.441 (average MW due to polydispersity of PEG)
Physical State	Powder
Storage	-20°C
Expiration Date	One year from date of receipt
M Lot Number	6313PJB011
Avanti Lot Number	880310P-25MG-B-011

ANALYSIS	SPECIFICATION	RESULTS
Physical examination	Powder: White solid which contains no foreign matter.	Pass
TLC (85:15 Chloroform:Methanol)	>99% Purity Ninhydrin: negative Iodine: one major spot Phosphorus: positive Charring: negative Water dip: one major spot	All Pass
HPLC-ELSD	NLT 99% (AUC)	100% (AUC)
FAME (GC/FID)	NLT 99% 14:0 (AUC)	100% (AUC)
Proton NMR	NMR spectrum consistent with structure	Consistent with structure

Approved By:



For Research Use Only. Not For Human Use.