

18:1 PEG3000 PE

1,2-dioleoyl-sn-glycero-3-phosphoethanolamine-N-[methoxy(polyethylene glycol)-3000] (ammonium salt)

880330C-25mg
Page 1 of 1

Molecular Weight	3770.622 (average MW due to polydispersity of PEG)
Physical State	Chloroform
Storage	-20°C
Expiration Date	One year from date of receipt
M Lot Number	6315CJA012
Avanti Lot Number	880330C-25MG-A-012

ANALYSIS	SPECIFICATION	RESULTS
Physical examination	Chloroform: Clear solution with no foreign matter.	Pass
TLC (85:15 Chloroform:Methanol)	>99% Purity Ninhydrin: negative Iodine: one major spot Phosphorus: positive Charring: positive Water dip: one major spot	All Pass

Approved By:

