

**18:1 PEG3000 PE**

1,2-dioleoyl-sn-glycero-3-phosphoethanolamine-N-[methoxy(polyethylene glycol)-3000] (ammonium salt)

**880330P-200mg**
**Page 1 of 1**

<b>Molecular Weight</b>	3770.622 (average MW due to polydispersity of PEG)
<b>Physical State</b>	Powder
<b>Storage</b>	-20°C
<b>Expiration Date</b>	One year from date of receipt
<b>M Lot Number</b>	6315PNB012
<b>Avanti Lot Number</b>	880330P-200MG-B-012

ANALYSIS	SPECIFICATION	RESULTS
Physical examination	Powder: White solid which contains no foreign matter.	Pass
TLC (85:15 Chloroform:Methanol)	>99% Purity Ninhydrin: negative Iodine: one major spot Phosphorus: positive Charring: positive Water dip: one major spot	All Pass

Approved By:

