

18:1 PEG3000 PE

1,2-dioleoyl-sn-glycero-3-phosphoethanolamine-N-[methoxy(polyethylene glycol)-3000] (ammonium salt)

880330P-25mg
Page 1 of 1

| | |
|--------------------------|--|
| Molecular Weight | 3770.622 (average MW due to polydispersity of PEG) |
| Physical State | Powder |
| Storage | -20°C |
| Expiration Date | One year from date of receipt |
| M Lot Number | 6315PJA012 |
| Avanti Lot Number | 880330P-25MG-A-012 |

| ANALYSIS | SPECIFICATION | RESULTS |
|------------------------------------|---|----------|
| Physical examination | Powder: White solid which contains no foreign matter. | Pass |
| TLC (85:15 Chloroform:Methanol) | >99% Purity Ninhydrin: negative Iodine: one major spot Phosphorus: positive Charring: positive Water dip: one major spot | All Pass |

Approved by:

