

14:0 PEG350 PE

1,2-dimyristoyl-sn-glycero-3-phosphoethanolamine-N-[methoxy(polyethylene glycol)-350] (ammonium salt)

880410C-25mg
Page 1 of 1

Molecular Weight	1019.287 (average MW due to polydispersity of PEG)
Physical State	Chloroform
Storage	-20 °C
Expiration Date	One year from date of receipt
M Lot Number	6324CJA012
Avanti Lot Number	880410C-25MG-A-012

ANALYSIS	SPECIFICATION	RESULTS
Physical examination	Chloroform: Clear solution with no foreign matter.	Pass
TLC (85:15 Chloroform: Methanol)	>99% Purity Ninhydrin: negative Iodine: one major spot Phosphorus: positive Charring: negative Water dip: one major spot	All Pass
Proton NMR	NMR spectrum consistent with structure	Consistent with structure

Approved By:

