

18:0 PEG1000 PE

1,2-distearoyl-sn-glycero-3-phosphoethanolamine-N-[methoxy(polyethylene glycol)-1000] (ammonium salt)

880720C-25mg
Page 1 of 1

Molecular Weight	1792.288 (average MW due to polydispersity of PEG)
Physical State	Chloroform
Storage	-20°C
Expiration Date	One year from date of receipt
M Lot Number	6340CJA025
Avanti Lot Number	880720C-25MG-A-025

ANALYSIS	SPECIFICATION	RESULTS
Physical Examination	Chloroform: Clear solution with no foreign matter	Pass
TLC (85:15 Chloroform: Methanol)	>99% Purity Ninhydrin: negative Iodine: one major spot Phosphorus: positive Charring: negative Water dip: one major spot	All Pass
FAME (GC/FID)	NLT 99% 18:0 (AUC)	99.8% (AUC)
HPLC-ELSD	NLT 99% (AUC)	99% (AUC)
Proton NMR	NMR spectrum consistent with structure	Consistent with structure

Approved By:



For Research Use Only. Not for Human Use.